

The SaHoMa®-II inhalation set for horses is equipped with the latest, modern technology:

> Highest acceptance due to the light construction and almost noiseless operation

> 275 x 250 x 325 mm, 1,200 grams total weight

Use of mesh technology for optimal particle sizes for treatment of the horses upper and lower airways

> 5,14 µm MMD | 3,92 µm MMAD approx. 45 % of the particles are  $\leq 5 \mu m$

Excellent aerosol delivery rate for a fast and effective inhalation therapy

0.8 ml NaCl / minute

## Technical features

The SaHoMa®-II inhalation set for horses is successfully used and recommended by many equine clinics, veterinarians, riding stables as well as professional and hobby equestrians:

## Horse clinics

Tierklinik Telate, Germany Tierklinik Lüsche in Barkum, Germany

### Showjumping

Paul Schockemöhle, Germany Ludger Beerbaum Stables, Germany Sportpferde Ehning, Germany

# **Eventing**

Michael Jung, German Olympic Squad Leonie Kuhlmann, German Junior Squad

#### Dressage

Hans Peter Minderhoud, Netherlands Gina Capellmann-Lütkemeier, Germany Dressurstall Rothenberger, Germany





Innovative Inhalationsgeräte

You would like to support your horse with an effective inhalation device?

We would be happy to assist you:



# **NEBU-TEC International GmbH**

Kreuzfeldring 17 D - 63820 Elsenfeld

Tel.: +49 (0) 6022 - 610 62 -0 Fax: +49 (0) 6022 - 610 62 -99 info@nebu-tec.de


shop.nebu-tec-shop.de

Follow us on

As a manufacturer of innovative and well-known inhalation devices for humans and animals, health is our top priority. We use our over 20 years of inhalation experience in human medicine to give your horse a healthy life. All the know-how and experience we have gained during this time flows into the development of our veterinary products and makes them unique.


www.nebu-tec.de

Instagram | Facebook | YouTube









Respiratory diseases are among the most common internal diseases in horses. The inhalation therapy with the SaHoMa®-II **inhalation set for horses** is the first choice for the treatment:

- Use for the treatment and prevention of inflammatory and chronic respiratory diseases in horses
- Support the respiratory system through humidification and care by using saline inhalation solutions
- **Inhaled application** of low-dose medications that **reach their** destination directly via respiration and develop their effect there better and faster than with systemic therapies - almost without any loss of active ingredients and undesired side effects
- The inhalation therapy reduces the need for visits to the vet and can be carried out by the horse owner on a long-term basis and without great effort
- Inhalation where, when and how your horse likes it best whether in the stable, during training on the lunge or during a leisurely walk








The field-tested inhalation set developed with veterinarians promises everything that is necessary for a quick and effective inhalation therapy. Compared to common inhalers, the difference of the SaHoMa®-II inhalation set for horses is in the details:





Integrated, automatic switch-off function after the filled inhalation solution has been completely nebulized

# Secure breathing


- Breathing without resistance for stress-free inhalation
- No contamination of the ambient air and no danger to surrounding living beings



**ADVANTAGES** 


- Wireless operation via a rechargeable, long-lasting battery integrated in the inhalation mask for unlimited mobility
- Stress-free inhalation for the owner thanks to an easy application and fast, hygienic cleaning
- Reliable application of medication and many other nebulization unit
- All-round feel-good factor through professional and individual support





Because we know that an inhaler is a big purchase for many horse owners, while they do not know how the horse reacts to the inhalation, we offer you to test the inhaler free of charge and without obligation for 14 days prior to purchase. Because only if you feel safe with the inhalation device and, above all, your horse can inhale comfortably and without stress, an effective inhalation therapy can take place. Interested? Please contact us for further information.





• Permanent monitoring of the inhalation process thanks to transparent construction Direct aerosol application in front

of the horse's nose for immediate inhalation with almost no loss of condensation and an efficiency of over 70 %

# Different sealing rings for fit and tightness

- Easy and inexpensive adaptation to all horse noses
- To seal the system and avoid losses